



FOCUS

Mexican Healthcare,
a System of Contrasts



AT A GLANCE

Key Considerations
for Entering China's
Healthcare Market



NON PROFIT

UNICEF on World AIDS Day

INDUSTRY NEWS



**MANUFACTURERS
LOOKING FOR
DISTRIBUTORS**

IEC INNOVATIVE
ENDOSCOPY
COMPONENTS, LLC
www.IECendoscopy.com



39-0105-00 -- Manual Camera Inspection System for
Endoscope and Instrument Repair, 1.0 mm & 2.0 mm
Camera probe available, 400 x 400 incl. 4 LED Lights,
L=3 meters, and 13" Flat panel screen

- Local and Webpage Video Storage
- Editing and Sharing Capabilities
- Multiple ways to organize your video library



Flexible Endoscope Inspection



Rigid Endoscope Inspection



Shaver System Inspection



Instrument Inspection

Endoscope CCD' Camera Repair Innovation

© INNOVATIVE ENDOSCOPY COMPONENTS, LLC



In an effort, to make a better and more functional product to serve the requests of our loyal customers, IEC has developed a new design for our small diameter endoscope CCD' repairs. This new design significantly reduces the diameter of the wire harness creating more room for other parts. It gives the endoscope CCD' the ability to angulate and move with less stress on the camera than ever before. We offer this service to endoscope repair companies and dealers worldwide.

Please contact us: Innovative Endoscopy Components, LLC: 320 International Parkway, Fort Lauderdale, FL 33325, USA

Tel: +1 954-217-8780,

Fax: +1 954-217-8781


► www.IECendoscopy.com

► info@endoscopy.md


→ **Smart Medical Fair**

www.smartmedicalfair.com/stand/h4p3b1z1553





VILLA



new
chronos
leading-edge R/F system

superior performance - operational flexibility

REDEFINED