

We can offer a worldwide visibility to your company.  
Even only one contract signed with one distributor of our database can help you multiply your ROI.

**3 2022**  
OCTOBER - DECEMBER



Edited by®  
Infodent  
Infomedix  
International  
Publishing &  
Consulting House



CHECK HERE WHO  
IS LOOKING FOR  
DISTRIBUTORS

VISIT US AT  
**MEDICA 2022**  
HALL 16, BOOTH A38



Worldwide  
Largest  
Inventory

New & Pre-Owned  
Endoscope  
Repair Parts



**AM350**

www.IECendoscopy.com  
320 International Pkwy  
Ft. Lauderdale, FL 33325  
954-217-8780



INFODENT Srl - Via dell'Industria 65 - 01100 Viterbo - Italy  
O/DCVT/h°5 fb del 24/05/02 - Sped. in A.P. - art. 1 comma 1 D.L. 353/2003 conv. In L. n. 46/04 - CDSU VT G.C.

## TOM GIANO HR **DC III**

### PERFECT.VISION



**DIRECT  
CONVERSION  
TECHNOLOGY:  
THE BEST X-RAY  
IMAGING QUALITY  
AVAILABLE TODAY,  
WITH ULTRA LOW DOSES.**

[www.newtom.it/en/medicale](http://www.newtom.it/en/medicale)

Cone Beam 3D Imaging  
**NewTom**  
what's next



IEC is the leading International ISO 9001/13485 certified endoscope, endoscope parts, equipment, specialty service provider, and consulting firm since 1999.



**Rigid Endoscopes**



**Endoscope Parts - Rigid**




**Endoscope Parts - Flexible**




**Injection Molding**




**Industrial Endoscopes**



**Tools & Devices**




**CCD' & Electronics Repair**



**Processor Panel Repair**



**Laser Marking Services**



**Contract Manufacturing**



**Training & Consulting**



**Stolen Endoscope Registry**

**Please contact us today for more information.**

320 International Parkway, Ft. Lauderdale, FL 33325 Tel (954) 217-8780 Fax (954) 217-8781 [www.IECendoscopy.com](http://www.IECendoscopy.com)



## Innovative Endoscopy Components, LLC

Premium Endoscope Repair Parts



**Innovative Endoscopy Components, LLC** has been the ISO13485/ISO9001 Certified vendor of choice to hundreds of endoscope service facilities and dealers worldwide, for over 24 years.

Our product range and services are constantly growing with international demand. Rapid prototyping, optical assemblies, injection molding, and CNC machined parts are offered just like OEM endoscope

and equipment labeling as well as CCD' repair and multilingual repair training and consulting.

Please contact us: Innovative Endoscopy Components, LLC: 320 International Parkway, Fort Lauderdale, FL 33325, USA Tel: 954-217-8780, Fax: 954-217-8781

[www.IECendoscopy.com](http://www.IECendoscopy.com)  
[info@endoscopy.md](mailto:info@endoscopy.md)



Visit us at: Smart Medical Fair 2022  
[www.smartmedicalfair.com/stand/h4p3b1z1553](http://www.smartmedicalfair.com/stand/h4p3b1z1553)



## Innovative Endoscopy Components, LLC

Endoscope CCD' Camera Repair Innovation



In an effort, to make a better and more functional product to serve the requests of our loyal customers, **IEC** has developed a new design for our small diameter endoscope CCD' repairs. This new design significantly reduces the diameter of the wire harness creating more room for other parts. It gives the endoscope CCD' the ability to angulate and move with less stress on the camera than ever before.

We offer this service to endoscope repair companies and dealers worldwide.

Please contact us: Innovative Endoscopy Components, LLC: 320 International Parkway, Fort Lauderdale, FL 33325, USA Tel: 954-217-8780, Fax: 954-217-8781

[www.IECendoscopy.com](http://www.IECendoscopy.com)  
[info@endoscopy.md](mailto:info@endoscopy.md)


Visit us at: Smart Medical Fair 2022  
[www.smartmedicalfair.com/stand/h4p3b1z1553](http://www.smartmedicalfair.com/stand/h4p3b1z1553)

## BMS Steril Maximum

Manufactured with the most up to date technologies, these are the main features:

- Stainless-steel chamber.
- Large graphic display, including a control software, making the autoclave extremely user-friendly.
- MTS (Memory Test System) a technical backup to assist the user with more reliable and faster diagnosis.
- Settings and menu in 19 different languages.
- Built-in printer, to carry out the daily data recording.
- Automatic electrical door lock for the safety of the user.
- 10 sterilization cycles: 6 Type-B Cycles and 4 Type-S Cycles according to the type of material to be sterilized.
- 2 test cycles ("Vacuum Test" and "Helix/Bowie&Dick Test"), see LINK.

**NEW!** The implementation of the conductivity meter for water control protects the steam generator, checking the quality of the water and preventing damage due to repeated use of bad quality water (this positively affect life cycle of the entire device).



**LINK:** to carry out the routine tests for the correct device functioning and sterilization efficiency, discover the quality of BMS dental sterilization control tests. Starting by the specific tests for Class B autoclaves, Bowie&Dick and Helix tests to the most used sterilization indicators: type 4 and type 6.

[www.bmsdental.it](http://www.bmsdental.it)  
[info@bmsdental.it](mailto:info@bmsdental.it)





Visit us at: Smart Medical Fair 2022, [www.smartmedicalfair.com/stand/h4p6b5z1522](http://www.smartmedicalfair.com/stand/h4p6b5z1522)



To be our partner it's easy  
just spin the wheel!

We are your partner  
for Endoscopes,  
Cameras,  
Shaver Blades,  
Repair Training  
& Consulting



We are specialized in the production  
of medical devices



**Silfradent srl**

Via Giuseppe di Vittorio 35/37 - 47018 S. Sofia (FC) - Italy  
+39 0543 970684 | [info@silfradent.com](mailto:info@silfradent.com) | [www.silfradent.com](http://www.silfradent.com)

**AM JSD**

[www.IECendoscopy.com](http://www.IECendoscopy.com)  
320 International Pkwy  
Ft. Lauderdale, FL 33325  
954-217-8780

